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en english

Dear patient,

Please read and follow the instructions for use carefully, and direct any questions to your physician, medical retailer, or our Customer Service Department at telephone number +49 (0)800 0 01 05 30 or info@bauerfeind.com.

Uses / Indications

Following individual adjustment by a professional, this orthotic blank is suitable for the conservative treatment of orthopedic foot and knee complaints and some hip and back pain.

Please note that your medical retailer will turn the orthotic blank into a foot orthosis that is individually tailored to your needs.

Orthopedic foot orthoses are used in the treatment of foot problems and foot deformities, such as flat foot, splay foot, pes valgus, mild cases of pes cavus, or a combination of these conditions. The treating physician is responsible for the precise diagnosis.

In addition to the treatment of foot complaints, foot orthoses can also be used to prevent excessive strain on the entire musculoskeletal system (knee, hip and back regions).

Risks of using this product

⚠ If there is not enough space inside the shoe for the foot and the selected foot orthosis, or if the product is worn in shoes with an incorrect length and width (or shoes with inadequate lateral guidance in the case of open shoe types), problems relating to poor fit may arise (e.g. a feeling of constriction or pressure, rubbing).

If the feet are positioned on the edges of the foot orthoses, blisters or calluses may form on the soles of the feet. The same symptoms may occur if excessive mechanical strain leads to holes or creases developing in the top covers, cushioning material or orthotic cores, or to the breakage of these. You must stop wearing foot orthoses if they are damaged in this way!

Creases in stockings or pantyhose must be pulled flat before the shoes are put on, as they can cause blisters and pressure sores.

If your symptoms increase, please contact your physician or medical retailer without delay.

The foot orthosis may not be used to treat patients with diabetic foot syndrome or rheumatoid foot deformities.

Application instructions

The ErgoPad work must be fitted to the feet and/or shoes by a professional¹.

AUSTRIA Bauerfeind Ges.m.b.H. Hainburger Straße 33 1030 Wien P +43 (0) 800 4430-130 F +43 (0) 800 4430-131 E info@bauerfeind.at	BENELUX Bauerfeind Benelux B.V. Waarderveldweg 1 2031 BK Haarlem THE NETHERLANDS P +31 (0) 23 531-9427 F +31 (0) 23 532-1970 E info@bauerfeind.nl	FRANCE Bauerfeind France S.A.R.L. B.P. 59258 95957 Roissy CDG Cedex P +33 (0) 1 4863-2896 F +33 (0) 1 4863-2963 E info@bauerfeind.fr	ITALY Bauerfeind Italia Srl Via Cornaggia 58 20092 Cinisello Balsamo (MI) P +39 02 8977 6310 F +39 02 8977 5900 E info@bauerfeind.it	SLOVENIA Bauerfeind d.o.o. Dolenjska cesta 242 b 1000 Ljubljana P +386 (0) 1 4272-941 F +386 (0) 1 4272-951 E info@bauerfeind.si	SPAIN Bauerfeind Ibérica S.A. Calle Albacete, 23, Entresuelo C 46007 Valencia P +34 96 385-6633 F +34 96 385-6699 E info@bauerfeind.es	SWITZERLAND Bauerfeind AG Vorderi Bode 5 5452 Oberrohrdorf P +41 (0) 56 485-8242 F +41 (0) 56 485-8259 E info@bauerfeind.ch	UNITED ARAB EMIRATES Bauerfeind Middle East FZ LLC Dubai Healthcare City Building 40, P.O. Box 505116 Dubai P +971 4 4335-684 F +971 4 4370-344 E info@bauerfeind.ae	USA Bauerfeind USA, Inc. 75 14 th St NE Suite 2350 Atlanta, GA 30309 P +1 800 423-3405 F +1 404 201-7800 E info@bauerfeindusa.com
BELGIUM P +32 (0) 2 527-4060 F +32 (0) 2 792-5345 E info@bauerfeind.be	BOSNIA AND HERZEGOVINA Bauerfeind d.o.o. Meše Selimovića 19 71000 Sarajevo P +387 (0) 33 710-100 F +387 (0) 33 619-422 E info@bauerfeind.ba	SERBIA Bauerfeind d.o.o. 102 Omladinskih brigada 11070 Novi Beograd P +381 (0) 11 2287-050 F +381 (0) 11 2287-052 E info@bauerfeind.rs	SINGAPORE Bauerfeind Singapore Pte Ltd. Blik 41 Cambridge Road #01-21 Singapore 210041 P +65 6396-3477 F +65 6295-3062 E info@bauerfeind.com.sg	CROATIA Bauerfeind d.o.o. Goleška 20 10020 Zagreb P +385 (0) 1 6542-855 F +385 (0) 1 6542-860 E info@bauerfeind.hr	GERMANY P +49 (0) 36628 66-3500 F +49 (0) 36628 66-3999 E info@bauerfeind.com BAUERFEIND.COM	ISRAEL P +972 (0) 3 523-1970 F +972 (0) 3 523-1970 E info@bauerfeind.co.il	NETHERLANDS P +31 (0) 23 531-9427 F +31 (0) 23 532-1970 E info@bauerfeind.nl	NETHERLANDS P +31 (0) 23 531-9427 F +31 (0) 23 532-1970 E info@bauerfeind.nl

Rev.17 – 2026-01_107544

When selecting the correct size of foot orthoses, the packaging information serves only as a guide. Just like correctly fitting shoes, foot orthoses must be at least 1.0 cm longer than the foot when the foot is in a weight-bearing, standing position **1**.

The shoes must guarantee sufficient internal volume for both the foot and the selected foot orthosis (adequate shoe length/width, removable, sufficiently thick cushioned sole, adequate lateral guidance of the foot orthoses in the case of open shoe types). An uncomfortable fit that constricts the foot must always be avoided **2**.

The forefoot (a1), midfoot (a2), and hindfoot (a3) must be well supported evenly from below by the foot orthosis and must not be positioned on one of its edges **3**.

Foot orthoses must always be worn as a pair! Any exceptions to this must be stipulated by your physician.

It must always be ensured that the left and right foot orthoses are clearly identified in a way which is easily understood by the patient.

The foot orthoses must not be exposed to temperatures exceeding 60 °C (e.g. heating, hair-dryer, direct sunlight).

The materials used for the foot orthoses are subject to aging and their properties can therefore change.

Cleaning instructions

Foot orthoses must be removed from the shoes on a daily basis in order to allow moisture to evaporate.

All surface materials facing the foot can be wiped if required, using a mild cleanser and a sponge or cloth. After cleaning, any residual soap must be thoroughly rinsed off with water. The foot orthoses cannot be washed in a washing machine.

The orthoses must be left to air-dry and no external heat source should be used (heating, hair-dryer, direct sunlight).

Mild disinfectants suitable for human skin may be applied sparingly (sprayed on). However, this may lead to the discoloration of top cover materials.

Permissible usage conditions / area of use

The foot orthoses must be worn in conjunction with hosiery in suitable shoes with a sufficient internal volume for both the foot and the selected foot orthosis.

Notes on reuse

Professionals will turn the orthotic blank into a medical device by modifying it as the indication requires. The medical device is designed to provide care for an individual patient and must not be used by any other person. In the event of improper modification or use, the medical effect will be compromised, and the guarantee as well as product liability restricted or excluded.

Foot orthoses are subjected to various types of mechanical strain (body weight, tension, pressure, torsion, friction, etc.) which affect the useful life of the foot orthoses. Temperatures, stocking materials, shoes, and environmental influences (e.g. perspiration) all play their part in reducing the useful life of the foot orthoses. After foot orthoses have been worn for a maximum of 6 months, we recommend their orthopedic functions be checked by a professional. The maximum wearing period of 6 months may reduce on a case-by-case basis, depending on the level of wear.

Warranty

The statutory regulations of the country of purchase apply. Where relevant, the country-specific warranty regulations between the distributor and buyer apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty.

Please do not make any alterations to the product yourself. This particularly applies to products that have been individually adjusted by a professional. Please follow our instructions for use and care. They are based on years of experience and ensure that our medically effective products continue to work for a long time. Only aids that function perfectly will help you. Failure to follow the instructions may also limit the warranty.

Duty to report

Due to regional legal regulations, you are obliged to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

You can dispose of the product in accordance with national legal requirements.

Modification instructions for fitting the ErgoPad work for safety shoes in accordance with DIN EN ISO 20345 and DIN EN 61340-5-1

The following materials must be used for further modification and bonding purposes:

Pryx® antistatic Purchase details: Fritz Minke GmbH & Co KG Cecilienstraße 31 47051 Duisburg Germany T: +49 (0) 203 281 01-0 F: +49 (0) 203 281 01-70 Email: gminke@minke.de	Renol = AL = (antistatic) Art. no. 24002601 (640 g) Purchase details: Ortho Planet Leder Brinkmann GmbH Eduard-Pestel-Strasse 3 49080 Osnabrück Germany T: +49 (0) 541 9 59 33-0 F: +49 (0) 541 9 59 33-33 Email: info@leder-brinkmann.de
---	---

Structural and modification specifications

The foot orthosis may only be modified using the procedure described below and in accordance with the defined product structure and modification instructions.

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