CONTACTS

AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

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Dear customer,
thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please read and observe these instructions for use carefully. If you have any questions, please contact your physician or medical retailer.

Intended purpose

CoxaTrain is a medical device. It is an orthosis¹ for joint stabilization and relief of pain in the hip joint for more movement with osteoarthritis of the hip.

Indications

- Osteoarthritis of the hip (Kellgren-Lawrence score 2–4)
- Muscular imbalance of the lumbar, pelvic and hip region
- Muscular problems after hip joint surgery
- Impingement of the hip joint
- Instability of the hip joint (after THR, after revision surgery)
- Feeling of instability in the hip
- Loosening of the hip prosthesis if surgery is not possible
- Chronic irritation (e.g. bursitis in conjunction with osteoarthritis)
- Chronic tendinopathy in conjunction with osteoarthritis

Risks of using this product

Please read and closely observe the specifications in these instructions for use and the instructions given by the specialist².

The benefits of this product are felt particularly during physical activity.

- Please remove your orthosis during longer periods of rest.
- Once you have been prescribed the CoxaTrain, only use it according to your indication and any additional instructions given by a medical specialist². Please consult a specialist or your physician if you want to use it together with other products. Do not modify the product yourself in any way; otherwise it may not help as expected or it may be harmful to your health. In these cases, any warranty or liability claims will be excluded.
- Avoid contact with ointments, lotions, or substances containing grease or acids.
- Any aids applied externally to the body may, if tightened excessively, lead to local pressure points or, in rare cases, constrict the underlying blood vessels or nerves. For this reason, check the fit of the product if you feel unpleasant pressure and consult the specialist in this regard or for general questions about handling.
- Caution: When loosening the straps or tensioning and fastening systems, or when removing the product or adjusting configured limitations of movement, there is not sufficient support / sufficient protection of the body part requiring treatment. In this case, particularly avoid placing any strain on this body part.
- If you notice any negative changes or an increase in symptoms while wearing the product, take the product off and contact your physician immediately.
- If you would like to use your Bauerfeind product due to acute symptoms / injuries, make sure to obtain and follow the advice of a medical professional before first using the product. Restrictions may exist given the specific indication which the specialist retailer can clarify. In particular, the ability to operate vehicles, other means of transport or machinery may be limited. In case of doubt, we recommend temporarily refraining from the activities mentioned above.

- Furthermore, proceed with extreme caution, as your ability to stand and walk is also heavily impacted.
- Check the specific ground characteristics and ensure the ground is not slippery whenever walking, standing or otherwise stressing the affected part of the body. If the ground is smooth, there is an acute risk of slipping and thus injury!

Contraindications

No clinically significant adverse reactions have been reported to date.

If you have any of the following conditions, the product should only be used after consultation with your physician:

- Skin disorders / injuries to the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up
- Impaired sensation, circulatory and movement disorders around the hips and legs
- Impaired lymph drainage – including soft tissue swellings of uncertain origin located away from the fitted aid

Components

- A – Hip cup with hinge
- B – Hip support
- C – Gluteal pad
- D – Flat trochanter cushion
- E – High trochanter cushion with hole
- F – Trochanter pad (white back)
- G – Trochanter pad (black back)
- H – Fabric joint protection
- I – Thigh support
- J – Covering cap
- K – 2 × limiting pins (available as an optional extra)
- L – SJ joint pads

Application instructions

Size determination (only by a trained specialist)

- The orthosis is selected and fitted by a specialist after measuring the patient. You can refer to the size chart on the packaging for the correct size depending on the measured values.
- To determine the right size, you will need to measure the pelvic circumference.
- You can determine the size of the thigh support based on the thigh circumference measured at the middle of the thigh. It is also necessary to indicate the affected side (left or right leg).
- To select the joint splint (short / long), you will need the patient's height.
s = short < 1.70 m height
l = long > 1.70 m height
In specific cases, it is possible that a "long" joint splint will fit better than a measured size of "short" or the other way round.

Adjustment of CoxaTrain (only by a trained specialist)

- First insert the hip cup with hinge  aligned according to the affected side into the welded pocket provided on the pelvic support for this purpose.  The Bauerfeind logo "B" on the hip cup should be facing outwards (away from the body) .
- Fasten the blue gluteal pad  (which is fastened on the right side by default) to the Velcro tab provided on the affected side (inner side of the pelvic support) .

- Then fasten the oval trochanter cushion  horizontally onto the triangular Velcro patches on the inner side of the hinge.
 - For targeted pressure relief of the greater trochanter, the second, perforated oval cushion  can be fastened at this point as an alternative .
 - After this, fasten the trochanter pad  (kidney-shaped, gray, white back) onto the joint splint above the trochanter cushion, keeping the hinge bent at a 90° angle .
 - The second trochanter pad included with delivery  (gray, kidney-shaped, black back) can be fastened below the pad  if required to generate an enhanced massage effect or to improve the action of the pad.
 - Now put the pelvic support with the assembled hinge on the patient. Make sure that the upper edge of the pelvic support is level with the upper edge of the iliac crest .
 - In case the hinge sticks out, you can reshape this to follow the anatomical progression of the thigh in the adjustable area of the splint  using your hands or the edge of the table up to 3 × .
 - To complete the product, first take off the support.
 - Then pull the fabric hinge protection  over the lower end of the hinge  and shorten using scissors if needed.
 - After this, insert the lower end of the hinge into the opening of the thigh knit  with the locking button pressed, then lock it in one of the six available positions (depending on the length of the leg) .
- TIP: Mark the number so the patient can find the right locking point after disassembly and reassembly during cleaning.
- Position the sacroiliac joint pads  (SI joint pads) vertically along the Velcro tabs so that the cut-out (hole) is situated at the posterior superior iliac spine (SIPS).
 - Depending on the indication / specification of the physician, the optionally available limiting pins  can be used to restrict bending and stretching in the hip joint (60°, 75°, 90°) .
 - To do so, position the limiting pin(s) in the corresponding holes (60°, 75°, 90°) .
 - Then snap the covering cap  onto the hinge .
 - Check the fit of the shaped orthosis on the patient and have the patient practice putting on and taking off the orthopedic aid in front of a mirror.
 - If the SI joint pads  get in the way of the pelvic support, you can unfasten them.

TIP: Make sure that the Velcro tabs do not touch the knitted fabric. This prevents damage to the material of the CoxaTrain.

Putting on the CoxaTrain (patient)

- Position the orthosis on the pelvis so that the upper edge of the knit fabric reaches no higher than the upper edge of the iliac crest with the hip cup and joint splint  aligned centrally on the side of the thigh .
- Then slide your fingers into the finger pockets provided on the fastening and steadily pull the fastening forward .
- Then press the left half of the fastening against the abdomen and position the right half of the fastening on top of the left half so that it can be fastened to it. When fastening the two halves,

remove your fingers from the left-hand finger pocket first and then from the right .

- Then take hold of the two tensioning straps and pull them forward, steadily and at the same time, until the desired level of tension (stabilization on the pelvis) is achieved .
- Then fasten the two tensioning straps on the large abdominal  fastening.
- Place the thigh knit  around the thigh by pulling it without creases  and close the Velcro fastening so that it sits flat . Alternatively, you can sit down while putting on the leg support for easier handling .

Taking off the CoxaTrain (patient)

- Open the Velcro fastening on the thigh knit .
 - Next, open the tensioning strap of the pelvic support   and attach it loosely to the Velcro tabs. Now you can open the fastening of the pelvic harness . To do so, hold the CoxaTrain firmly on the side with the hip cup and hinge .
- TIP: For easier handling and to prevent damage to the material of the CoxaTrain, always fasten the opened Velcro fastenings directly to the velour fabric of the orthosis.

Cleaning instructions (patient)

- First please note the side on which the locking button is locked into the leg knit. There are six available positions (depending on the length of the leg) .
- When adjusting the orthosis, a mark should be placed for better orientation.
- Remove the thigh knit from the joint splint by pressing the locking button down and pulling the thigh knit out of the joint splint .
- Please wash the orthosis, pads and cushion by hand with a mild detergent at 30°C.
- Wipe the hip cup and aluminum parts with a soft cloth and a mild cleaning agent.
- Then leave the orthosis to air dry.
- After drying, reattach the thigh knit  to the joint splint  by pressing the locking button and sliding the thigh knit with the thigh plate onto the joint splint .
- Lock the knit into its previous position .
- Fasten the pads back into the pelvic support or on the affected side as pictured .
- Regular care guarantees optimum effectiveness of the product. Also observe the instructions on the label.
- Never expose the product to direct heat as this may impair its effectiveness.
- This product has been tested by our own quality control system. Nevertheless, if you have a complaint, please contact your specialist on site.

Part of the body this product is used for Hip joint

Assembly and fitting instructions

The product is to be adjusted by specialists.

Maintenance instructions

If handled and cared for correctly, the product is maintenance-free.

Technical specifications / parameters

Your product comprises:

- Knitted mesh
- Velcro fastenings
- A lateral plastic splint with aluminum hinge
- Pads / cushions

- Limiting pins (available as an optional extra)

Notes on reuse

This product is intended to be fitted to one patient only.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the CoxaTrain have not been observed.

Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the specialist were not observed
- The product was modified arbitrarily.

Disclaimer

Do not diagnose yourself or decide on which medication to take unless you are a medical specialist. Before first using our medical product, please seek the advice of a physician or trained specialist, as this is the only way to assess the effect of our product on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of this specialist as well as all instructions in this document or its online version, including any excerpts (as well as texts, images, graphics etc.). If you still have doubts after your consultation with the specialist, please contact your physician or retailer, or contact us directly.

Duty to report

Due to regional legal regulations, you are obliged to immediately report any serious incident involving the use of this medical product to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Pelvic support with SI joint and gluteal pads

Polyamide (PA), Styrene-ethylene / butylene-styrene block copolymer (TPS-SEBS), Polyurethane (PUR), Polypropylene (PP), Elastane (EL), Cotton (CO), Polyester (PES), Polyoxymethylene (POM), Ethylene-vinyl acetate (EVA), Polyethylene (PE), polyacrylate adhesive

Thigh support

Polyamide (PA), Polyurethane (PUR), Elastane (EL), Glass fiber (GF), Polyester (PES)

Hinge with trochanter pad + cushion

Polyamide (PA), Aluminum (Al), Glass fiber (GF), Polyketone (PK), zinc alloy, Polyoxymethylene (POM), Styrene-ethylene / butylene-styrene block copolymer (TPS-SEBS), steel, galvanized, Cotton (CO), Polyurethane (PUR), Ethylene-vinyl acetate (EVA), steel, stainless, Polyethylene (PE), Polyester (PES), Elastane (EL), highly abrasion-resistant plastic, polyacrylate adhesive, Polytetrafluorethylene (PTFE)

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s = short <1,70 m de taille corporelle

l = long >1,70 m de taille corporelle

¹ Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso

² A specialist is any person who is authorized according to the state or institutional regulations as well as other provisions relevant to them for fitting and instruction in the use of orthopedic aids. As the manufacturer, we strongly recommend that this person be trained in orthopedics or possess comparable acquired competencies as well as dispose of orthopedic infrastructure.

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