
BAUERFEIND AG
 Triebeser Straße 16
 07937 Zeulenroda-Triebes
 Germany
 P +49 (0) 36628 66-40 00
 F +49 (0) 36628 66-44 99
 E info@bauerfeind.com
BAUERFEIND.COM

RhizoLoc® OA

Daumenorthese | Thumb orthosis | Orthèse de pouce
 Ortesi per il pollice | Ortesis de pulgar

Gebrauchsanweisung | Instructions for use
 Notice d'utilisation | Istruzioni per l'uso
 Instrucciones de uso

1/2
2/2

MADE
 IN
 GERMANY

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(en) english

Dear Customer,

thank you for choosing a Bauerfeind product. We work to improve the medical effectiveness of our pro-

ducts every day – because your health is very important to us. Please read and observe these **instructions for use** carefully. If you have any questions, please contact your physician or medical retailer.

Intended purpose

The RhizoLoc OA is a medical device¹. It is an orthosis¹ for the stabilization and relief of the carpometacarpal joint and metacarpophalangeal joint of the thumb.

Indications

- Osteoarthritis (e.g. osteoarthritis of first carpometacarpal joint)
- Arthritis
- Acute / post-traumatic (e.g. excessive strain, sprains, contusions)
- Chronic (ligament instability, recurring joint misalignments)
- Post-operative (e.g. after arthroplasty)
- Functional conditions (excessive strain, in case of recurring joint misalignments etc.)
- Irritation (chronic, acute, post-traumatic, post-operative)

Risks of using this product

⚠ Caution*

Please closely observe the specifications in these instructions for use and the instructions given by the specialist.

- RhizoLoc OA must only be worn in accordance with the specifications contained in these instructions and used for the areas of application (places of use) listed.

- Correct usage / fitting is assumed.
- To guarantee the best possible fit of the RhizoLoc OA, it is essential to check whether the product is adequately shaped and whether it needs to be anatomically reshaped by a specialist². Only a trained specialist may carry out the (initial) fitting of the orthosis and provide instruction.
- The RhizoLoc OA label with information about the size, side specification and CE marking can be found on the product packaging and on the inside of the thumb brace **(A)**.
- Use in combination with other products, e.g. as part of compression treatment (arm stockings), must be discussed beforehand with your physician.
- No product liability is accepted in the event of improper use, and the legal guarantee is limited or voided entirely.
- Do not allow the product to come into contact with ointments, lotions, or substances containing grease or acids. Do not make any alterations to the product.
- Side effects involving the body as a whole have not been reported to date.
- Do not put the RhizoLoc OA too tightly as this could lead to local pressure symptoms. In rare cases, constriction of blood vessels and nerves may occur. If this happens, loosen the straps of the RhizoLoc OA and, if necessary, check the size of your RhizoLoc OA.
- If you notice any negative changes or an increase in symptoms while using the

product, stop any further use and contact your physician.

- **Caution:** When loosening the straps or tensioning and fastening systems, or when removing the product or adjusting configured limitations of movement, there is not sufficient support / sufficient protection of the body part requiring treatment. In this case, particularly avoid placing any strain on this body part.
- If you would like to use your Bauerfeind product due to acute symptoms / injuries, make sure to obtain and follow the advice of a medical professional before first using the product. Restrictions may exist given the specific indication which the specialist retailer can clarify. In particular, the ability to operate vehicles, other means of transport or machinery may be limited. In case of doubt, we recommend temporarily refraining from the activities mentioned above.

Contraindications

No clinically significant hypersensitive reactions have been reported to date. If you have any of the following conditions, such aids should only be fitted and worn after consultation with your physician:

- Skin disorders / skin lesions in the relevant part of the body, particularly if inflammation is present; also in case of any raised scars with swelling, redness or excessive heat build-up
- Impaired sensation and circulatory disorders in the

arm / hand (e.g. diabetes mellitus)

- Impaired lymphatic drainage – including unspecified soft tissue swelling located in distal direction of the fitted orthosis
- Mental and physical disabilities that render it impossible to use the orthosis safely

Components

- A** – Thumb brace
- B** – Long strap (wrist)
- C** – Self-adhesive (textile) wraparound thumb strap, a spare strap is enclosed
- D** – Fastening with eyelet
- E** – Strap buckle with prong
- F** – Short strap (palmar)
- G** – Strap connector

Application instructions

Before using the orthosis for the first time, practice taking it on and off with a medical specialist.

Putting on the RhizoLoc OA

To put the RhizoLoc OA on the affected thumb, the fastenings **(B)** and **(C)** are opened. Guide the hand through the wrist strap **(B)** with the thumb extended to the side until the plastic brace is positioned between the thumb and index finger. **1** In the process, you can also use the hand symbol on the inside of the thumb brace for orientation **(A)**.

Now guide the long strap **(B)** around your wrist. The fastening with eyelet **(D)** hooks into the strap buckle on the back of the hand **(E)**. You will hear a soft clicking sound. **2** Now wrap the self-adhesive wraparound thumb strap **(C)** around your thumb without pulling tightly. **3**

Note: If the textile wraparound thumb strap **(C)** becomes wet, it should be removed and allowed to air-dry. Please use the enclosed spare strap while the removed strap is drying. To do so, unclick the strap with rivet directly out of the orthosis **4** and engage the enclosed second strap at the same position on the inside of the orthosis. Once the wet strap has dried, it is ready to be used again.

Following consultation with your physician, you can loosen the stabilization of the thumb basal joint depending on the phase of your treatment by taking off the wraparound thumb strap **4**. This enables a moderate dynamization of the thumb basal joint.

Taking off the RhizoLoc OA

First unhook the long strap **(B)** with the fastening **(D)** from the prong of the strap buckle **(E)** and then unwrap the wraparound thumb strap **(C)** – again without pulling tight. Now take off the RhizoLoc OA. **1**

Cleaning instructions

Note: Never expose the product to direct heat (e.g. heaters, sunlight, in your car, etc.). In case of contact with salt water or chlorinated water, rinse the orthosis out with plain water. This may cause damage to the material. Which can impair the effectiveness of the product.

1. Please close all fastenings and only wash the orthosis¹ by hand at 30°C using a mild detergent or cleaning fluid. All components can stay on the orthosis during washing.

2. The self-adhesive wrap-around thumb strap (**D**) will not adhere when moist, if this is the case swap the strap with the enclosed second wraparound thumb strap. Once the wet strap has dried, it can be used again.
3. Regular care guarantees optimum effectiveness of the product.
4. Do not tumble-dry the product. Leave the orthosis to air dry.

Part of the body this product is used for

In accordance with the indications (Thumb). Please see intended purpose.

Maintenance instructions

If handled and cared for correctly, the product is maintenance-free.

Assembly and fitting instructions

The product is to be adjusted by specialists².

Technical specifications / parameters

Your product comprises:

- Thumb brace (**A**) with a long (wrist) strap (**B**) and a short (palmar) strap (**F**)
- Wraparound thumb strap (**C**) (optional)
- Strap connector (**G**)
- Fastening with eyelet (**D**) and strap buckle (**E**)

Notes on reuse

This product is intended to be fitted to one patient only.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the RhizoLoc OA have not been observed.

Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the specialist were not observed
- The product was modified arbitrarily.

Disclaimer

Do not diagnose yourself or decide on which medication to take unless you are a medical specialist. Before first using our medical product, please seek the advice of a physician or trained specialist, as this is the only way to assess the effect of our product on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of this specialist as well as all instructions in this document or its online version, including any excerpts (as well as texts, images, graphics etc.). If you still have doubts after your consultation with the specialist, please contact your physician or retailer, or contact us directly.

Duty to report

Due to regional legal regulations, you are obliged

to immediately report any serious incident involving the use of this medical product to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Thermoplastic elastomer (TPE), Polyethylene, low-density (PE-LD), Polyester (PES), Thermoplastic polyurethane (TPU)

 – Medical Device

 – DataMatrix Unique Device Identifier UDI

Specialists²

Fitting of strap connector (G)

For immobilization and relief of the carpometacarpal joint and metacarpophalangeal joint of the thumb, the thumb brace (**A**) is positioned between the thumb and index finger.  The long wrist strap (**B**) is guided over the back of the hand.  The short strap (**F**) lies on the inside of the hand.

Guide the long wrist strap around the wrist and check that the orthosis and straps fit closely. **Note:** During fitting, make sure that the hand and thumb are in a functional position (see cover sheet of instructions for use).

If the straps are positioned too loosely/too firmly or not closely fitted to the hand, this can be adjusted.

To do so, take off the orthosis and remove the fastening with eyelet (**D**) and the strap connector (**G**) completely from the long wrist strap.

Reposition the thumb brace between the thumb and index finger  and . Guide the long wrist strap again from the back of the hand around the wrist so it lies flat. At the crossing point between the short strap (**F**) and the long strap , mark the two overlapping holes.

 Now position the strap connector (**G**) in the marked hole of the short (palmar) strap (**F**). Note: The smooth surface of the strap connector should be facing the hand.

 Next insert the long wrist strap (**B**) back into the strap connector (**G**) and pull the strap until the two marked holes overlap.

TIP:  In case the short strap (**F**) protrudes over the long strap (**B**), the endpiece of the short strap (**F**) can be shortened with scissors.

Adjustment of strap buckle with prong (E)

 By default, the strap buckle (**E**) is already fastened to the wrist strap (**B**). Ideally it should be shifted to the height of the index and middle finger.

Once the bracket is positioned properly between the thumb and index finger, next guide the long wrist strap to the back of the hand again and check whether the strap passes over the prong of the strap buckle. If this is not the case, the strap buckle (**E**) can be moved. There are seven possible positions in total. The ideal fastening

position (with the long strap in place) is when the orthosis lies flat, is comfortable and is not too tight on the hand / on the wrist.

Length of the wrist strap (B)

① Now hook the long strap without the fastening (D) into the strap buckle. If the orthosis lies flat, is comfortable and is not too tight on the hand, the fastening position is ideal. Now the strap needs to be shortened. To do so, cut the hooked hole as well as the hole beneath it off the strap in a semicircular shape.

② Click the fastening with eyelet (D) into the last two holes of the long wrist strap (B).

Please note the **instructions for use** for patients.

Note: If the self-adhesive wraparound thumb strap (C) is too long, this can also be shortened using scissors.

Version: 2024-03

¹ Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso

² A specialist is any person who is authorized, according to their state's regulations, to fit and provide instructions on the use of supports and orthoses.

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Bauerfeind France S.A.R.L.