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Dear patient,
Thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please closely observe the specifications in these instructions for use and the instructions given by the professional².
If you have any questions, please contact your physician or medical retailer.

Intended purpose
SecuTec Genu Flex is a medical device. It is a stabilizing orthosis¹ that limits the range of motion in the case of complex knee injuries.

- Indications**
- Rupture of the anterior and/or posterior cruciate ligament (ACL / PCL)
 - After ligament surgery / ligament reconstruction
 - Severe and/or complex instability (traumatic, degenerative, e.g. "unhappy triad")
 - Collateral ligament injuries
 - To restrict the knee joint's range of motion (e.g. after meniscus refixation or meniscal implantation)
 - For non-surgical treatment or long-term rehabilitation
 - Functional instabilities

- Risks of using this product**
- Only use the SecuTec Genu Flex according to your indication and the specific instructions given by a professional².
 - Before using the product, the professional² or your physician must determine whether thrombosis prevention is necessary.
 - Please consult a professional² or your physician if you want to use it together with other products. Do not modify the product yourself in any way; otherwise it may not help as expected or it may be harmful to your health. In these cases, any warranty or liability claims will be excluded.
 - Avoid contact with substances containing grease or acids, ointments, lotions, extreme temperatures as well as chemicals with an oxidizing effect and strong acids.
 - Avoid contact with abrasive or chafing substances (such as sand) to maintain the functionality of the orthosis.
 - Any aids applied externally to the body may, if tightened excessively, lead to local pressure points or, in rare cases, constrict the underlying blood vessels or nerves.
 - If you notice any changes or an increase in symptoms while wearing the product, stop any further use and contact your physician.
 - Caution: When loosening the straps, or when doffing the product or adjusting configured limitations of movement, there is not sufficient protection of the body part requiring treatment. In this case, particularly avoid placing any strain on this body part.
 - If you become aware of any changes to the product while using it, consult a professional².
 - If you want to use the SecuTec Genu Flex more than 24 months upon receipt, we recommend having a professional² check the plastic of the orthosis for damage after this period because they can prematurely age due to sunlight, heat, cold and moisture, for example. If your usage habits or intensity of usage differ from the therapy instructions given by your treating physician or from the therapy guidelines, the product's period of use may be limited.
 - Should any problems with the fit of the orthosis arise, for instance if the shape of your leg changes (e.g. due to an increase or reduction in muscle), or should the product have a functional defect, please immediately contact your physician or professional².
 - Side effects involving the body as a whole have not been reported to date. Correct usage / donning is assumed.
 - Please adjust your behavior to the situation associated with your specific condition and any associated restriction of movement, and avoid anything that could delay recovery where possible. Proceed with particular caution.
 - Furthermore, proceed with extreme caution, as your ability to stand and walk may be impacted.
 - Ask your physician or professional² which specific restrictions are personally relevant for you and must be observed.
 - If your perceptive faculties are restricted, if you take painkillers or blood-thinning medications or suffer from dementia, our products may only be worn when prescribed by a physician.
 - If you would like to use your Bauerfeind product due to acute symptoms / injuries, make sure to obtain and follow the advice of a professional before first using the product. Restrictions may exist given the specific indication which the professional² can clarify. In particular, the ability to operate vehicles, other means of transport or machinery may be limited. In case of doubt, we recommend temporarily refraining from the activities mentioned above.

- Contraindications**
- No clinically significant adverse reactions have been reported to date. If you have any of the following conditions, the product should only be used after consultation with your physician:
- Skin disorders / injuries to the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up
 - Varicose veins (varicosis)
 - Impaired sensation and circulatory disorders of the legs / feet (e.g. diabetes mellitus)
 - Impaired lymphatic drainage – including unspecified soft tissue swelling located in distal direction of the fitted aid

Application instructions
Size determination:
The size is selected and fitted by a professional² as needed after measuring the patient. You can refer to the size chart on the packaging for the correct size.

Procedure for the first fitting of the SecuTec Genu Flex knee orthosis
The connecting element with the B logo defines the front side of the orthosis. It must be positioned on the front of the thigh. In addition, the upper end of the orthosis is marked by arrows on the frame. On Bauerfeind products, the arrows always point upwards.

- Adjustment options for the hinge (only by a professional²)
- The hinge cover can be opened via a screw fastening. This cover fastens the movement-limiting wedges (extension – front / flexion – rear) and also protects the hinge mechanism against dirt.
 - Opening the screw fastening (quarter turn until the arrow points to ) allows the hinge protection to be removed and the desired limiting wedges to be changed.
 - The desired extension and flexion wedges are selected from the enclosed frame and inserted  .
 - The following limitations are possible:
Extension: 0°, 10°, 20°, 30°, 45°
Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
 - Then, click the hinge cover onto the hinge, check whether it is correctly positioned and firmly attached to the hinge, and close the screw fastening  . Make the adjustment to both hinges.

- Adjustment of the frame (only by a professional²)
- To ensure an optimal fit, the hinges must be positioned parallel to the shin bone and the thigh. The middle joining elements are positioned centrally on the calf and the thigh. .
 - The circumference around the proximal and distal frame of the SecuTec Genu Flex can be adjusted by removing and adding the Flex-Links.
 - To do this, open the frame closure with a pointed tool .
 - Separate the Flex-Links from the pad by pulling carefully on the pad and thereby removing the plastic snap fasteners from the frame  a–b.
 - Now you can loosen the connection between the Flex-Links by bending against the curvature of the frame .
 - To do so, grab the two Flex-Links that you want to separate and slightly bend them against the curvature of the frame to loosen the outer two guide edges. Then bend the Flex-Links by 90° to each other to loosen the center piece. You will hear a clicking sound.
 - As needed, you can now either remove the Flex-Links or add new parts from the packaging included.
 - Each chain link corresponds to a circumference change of 1 cm.
 - To join together the Flex-Links, first slide the two outer guide edges into the groove of the front Flex-Link  a.
 - After this, the center piece can be clicked in by bending slightly against the curvature of the frame  b.
 - It is recommended to carry out alteration on the inside and outside of the orthosis symmetrically to the middle joining element.
 - To reconnect the frame with the Flex-Modules, slide the frame closure into the frame mount .
 - Finally, fasten the pad back on the frame using the snap fasteners  a–b.
- Proceed symmetrically starting from the middle joining element.

Adjustment of the orthosis height (only by a professional²)
You can adjust SecuTec Genu Flex either distally or proximally in two lengths.
To do so, press a pointed tool against the plastic button in the area of the frame mount  and slide the hinge leg upwards until it engages in the next position. Standard delivery is made in the shortest setting.

Crimping the hinge legs (only by a professional²)
If required, you can reshape the hinge leg using round crimping irons with protective caps. When doing so, the plastic parts can remain on the hinge legs.

Static fitting (only by a professional²)

- Perform the initial adjustment with the straps open.
- Bend the knee to be treated to a 30° angle and place the SecuTec Genu Flex on it so that the hinge pivot point of the orthosis hinge is aligned with the compromise axis of rotation as defined by Nietert .

- Donning of the SecuTec Genu Flex**
- Don the orthosis against your skin. Wearing it over pants can cause slipping and impair functionality.
 - The leg is bent at a gentle angle (approximately 30°). The position of the kneecap is then determined through palpation.
 - The center of the hinge (pivot point of the orthosis hinges) is positioned at approximately the same height as the center of the kneecap .
 - When tightening the straps, follow the steps from I to IV . Start by pulling the strap tight directly beneath the knee, then proceed to the straps closer to the joint above the knee. Then tighten the straps further from the joint in the same order. When closing the straps, check whether the orthosis has rotated laterally around the leg and correct this if necessary.
 - After all the straps have been closed, the orthosis must exert an even pressure when the muscles are tensed.
 - Once the orthosis has been put on, check once more that it is in the correct position.

Part of the body this product is used for In accordance with the indications (knee). Please see intended purpose. Assembly and fitting instructions The SecuTec Genu Flex is supplied as a pre-assembled orthosis in standard sizes. The first fitting, which includes the hinge adjustments along with the adjustment of the frame and orthosis height, crimping and static fitting, must be carried out by a professional². Maintenance instructions If handled and cared for correctly, the product is maintenance-free. TecÜical specifications / parameters The SecuTec Genu Flex is an orthosis for the knee joint. It consists of a frame with individual Flex-Links, hinges, straps and cushions. Notes on reuse This product is intended to be fitted to one patient only.	Warranty The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the SecuTec Genu Flex have not been observed. Warranty is excluded if: • The product was not used according to the indication • The instructions given by the professional² were not observed • The product was modified arbitrarily Notice of liability For use of the SecuTec Genu Flex, the following applies: Do not diagnose yourself or decide on which medication to take unless you are a professional².
---	---

- Doffing of the SecuTec Genu Flex**
- When doffing the orthosis, it is recommended to adopt a seated position.
 - Open the straps marked in the image  in reverse order (IV to I) and slide the orthosis down on the leg. Do not open the remaining two straps to ensure a uniform fit. TIP: If you have been strictly instructed not to bend your leg, ask someone around you to help pull the orthosis over your foot.

Cleaning instructions
Never expose the orthosis to direct heat (e.g. heater, sunlight, etc.)!
• The aluminum frames of the orthosis are coated and can be cleaned using pH-neutral soap like the plastic parts.
• We recommend washing the pads and straps by hand at 30 °C. To reinsert the straps, proceed as shown in images  a–d.
• After contact with salt water and chlorinated water, thoroughly rinse the orthosis out with plain water.
• If you notice that a hinge is stiff after contact with coarse dirt (e.g. with sand), contact your professional² to have the hinge cleaned. This product has been tested by our own quality control system. However, should you have any complaints, please contact your local specialist orthopedic store or medical supply retailer.

Before first using our product, please seek the advice of this professional², as this is the only way to assess the effect of our SecuTec Genu Flex on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of the professional² as well as all instructions in this brochure or its online version, including any excerpts (as well as texts, images, graphics etc.). If you still have doubts after your consultation with the professional², please contact your physician or distributor, or contact us directly. Duty to report Due to regional legal regulations, you are obliged to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.	Disposal Realignment orthosis: Once you have finished using the product, please dispose of it according to local specifications, send the product back using the included return slip or return the product to the medical retailer who supplied you with the orthosis. Material content Polyketone (PK), Aluminum (Al), Polyamide (PA), Polyurethane (PUR), steel, polished, anodized, Polyester (PES), Polyoxymethylene (POM), Thermoplastic polyurethane (TPU), Polypropylene (PP), Elastane (EL), stainless steel, Glass fiber (GF), Polyethylene terephthalate (PET), brass, nickel-free, adhesive gum resin  – Medical device  – DataMatrix Unique Device Identifier UDI
---	---

¹ Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso.
² A professional is any person who is authorized according to the state or institutional regulations as well as other provisions relevant to them for fitting and instruction in the use of orthopedic aids. As the manufacturer, we strongly recommend that this person be trained in orthopedics or possess comparable acquired competencies as well as dispose of orthopedic infrastructure.

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