AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

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 english

Dear patient,

Thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please read and observe these instructions for use carefully. If you have any questions, please contact your physician or medical retailer.

Intended purpose

SecuTec Lumbo is a medical device. It is an orthosis¹ for lordosis reduction and relief of the lumbar spine with a bridging frame.

Indications

- Most severe lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Extremely severe (pseudo) radicular lumbar syndrome
- Spondylolysis / spondylolisthesis, grades I–III
- Moderate to severe lumbar deformity in cases of facet syndrome / osteoarthritis
- Lumbar spinal canal stenosis
- Moderate to severe degeneration / muscular insufficiency
- Vertebral fractures (lumbar spine) with vertebral body damage
- Condition following herniation, e.g. post-operative
- Condition following discectomy

Risks of using this product

Caution³

- SecuTec Lumbo must only be worn in accordance with the specifications contained in these instructions for use and used for the areas of application (places of use) listed.
- SecuTec Lumbo is a prescribed product that should be applied under a physician's guidance.
- To ensure the best possible fit, the SecuTec Lumbo must be individually adapted by a professional² to optimally suit the patient's requirements (adaptation of the bridging frame).
- Only a professional² may carry out the (initial) fitting of the product and provide instruction.
- Use in combination with other products, e.g. as part of compression treatment (compression stockings), must be discussed beforehand with your physician.
- The SecuTec Lumbo is to be worn over your T-shirt or undershirt.
- Side effects involving the body as a whole have not been reported to date. Correct usage / donning of the product is assumed. Do not don the SecuTec Lumbo too tightly as this could lead to local pressure symptoms. In rare cases, constriction of blood vessels and nerves may occur.
- If you notice any negative changes or an increase in symptoms while using the product, stop any further use and contact your physician.
- No product liability is accepted in the event of improper use.
- Caution: When loosening the straps or tensioning and fastening systems, or when doffing the product or adjusting configured limitations of movement, there is not sufficient support / sufficient protection of the body part requiring treatment. In this case, particularly avoid placing any strain on this body part.
- Do not allow the product to come into contact with ointments, lotions, or substances containing grease or acids.
- If you would like to use your Bauerfeind product due to acute symptoms / injuries, make sure to obtain and follow the advice of a professional before first using the product. Restrictions may exist given the specific indication which the professional can clarify. In

particular, the ability to operate vehicles, other means of transport or machinery may be limited. In case of doubt, we recommend temporarily refraining from the activities mentioned above.

 Contains natural rubber latex (elastodiene), which can cause skin irritation in persons with an allergy. Doff the orthosis if you notice any allergic reaction.

Contraindications

No clinically significant hypersensitivity reactions have been reported to date. If you have any of the following conditions, such aids should only be donned and worn after consultation with your physician:

- Skin disorders / injuries to the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up
- Impaired sensation and movement disorders, e.g. in the case of diabetes mellitus
- Severe restriction of the performance of the heart and lungs (risk of a rise in blood pressure if the aid is worn during intense physical exertion)

Application instructions

A – Knitted fabric portion
B – Fastening with finger pocket
C – Tensioning straps
D – Pocket with bridging frame
E – Abdominal pad

Donning of the product

The SecuTec Lumbo must be adapted by a professional² before the first donning.

❶ Position the orthosis centrally, at waist height, so that the lower back is covered. To ensure optimal effectiveness, the orthosis should be positioned as low as possible without pressing uncomfortably against the groin area while the patient is sitting upright. **❷** Then slide your fingers from the sides into the finger pockets provided on the fastening and steadily pull the fastening forward.

❸ Then press the left half of the fastening onto your stomach and position the right half of the fastening on top of the left half so that it can be fastened to it on the left-hand side. When fastening the two halves, remove your fingers from the left-hand finger pocket first and then from the right. **❹** Then take hold of the two tensioning straps and pull them forward, steadily and at the same time, until the desired level of tension (stabilization) is achieved. **❺** Then fasten the two tensioning straps close together on the large fastening. **❻ + ❼** Now fasten the two upper tensioning straps.

Doffing of the product

To doff the orthosis, start by undoing the four tensioning straps. Loosely refasten the straps onto the two halves of the fastening. Then open the large fastening. Doff the orthosis, then refasten the fastenings to avoid damage.

Cleaning instructions

Note: Never expose the SecuTec Lumbo to direct heat (e.g. heater, sunlight, in your car etc.)! This may cause damage to the material, which can impair the effectiveness of the SecuTec Lumbo.

The sewn-in label of the SecuTec Lumbo, containing information about the product name, size, manufacturer, care instructions and CE identification, is located on the inside of the knitted material.

Before washing, remove the bridging frame from the velour pocket, unfasten the abdominal pad and close all Velcro fastenings to avoid damage.

Wipe the bridging frame clean with a mild detergent. Wash the knitted fabric portion and the abdominal pad separately by hand with a mild detergent at 30 °C.

Allow the bridging frame, abdominal pad and knitted section to air dry.

When reinserting the bridging frame into the velour pocket provided, make sure that the two markings “up” and outward are correctly oriented.

Part of the body this product is used for

In accordance with the indications (lower back). Please see intended purpose.

Maintenance instructions

If handled and cared for correctly, the product is practically maintenance-free.

Assembly and fitting instructions

The SecuTec Lumbo is supplied in a standard size and ready for use. The bridging frame must be individually adapted by a professional and positioned in the knitted section.

Technical specifications / parameters

The SecuTec Lumbo is an orthosis for the lower back. It consists of a knitted section, a velour pocket for inserting the bridging frame, an abdominal pad and six double-sided micro Velcro fastenings. The velour pocket and bridging frame are available as spare parts on request.

Notes on reuse

The SecuTec Lumbo is intended for your treatment only. It has been adjusted specifically for you. For this reason, do not pass the SecuTec Lumbo on to other people.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. The warranty may be limited or excluded if the instructions on how to handle and care for the SecuTec Lumbo have not been observed. Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the professional were not observed
- The product was modified arbitrarily

Disclaimer

Do not diagnose yourself or decide on which medication to take unless you are a professional². Before first using our medical device, please seek the advice of a physician or professional² as this is the only way to assess the effect of our product on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of this professional as well as all instructions in this document or its online version, including any excerpts (as well as texts, images, graphics etc.). If you still have doubts after your consultation with the professional, please contact your physician or distributor, or contact us directly.

Duty to report

Due to regional legal regulations, you are obliged to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Polyamide (PA), Polyurethane (PUR), stainless steel, Viscose (CV), Glass fiber (GF), Elastodiene (ED), Polyethylene (PE), Polyoxymethylene (POM), polyether foam, Ethylene propylene diene monomer rubber (EPDM), Polyester (PES), Carbon fiber (CF), Cotton (CO), Polypropylene (PP), Elastane (EL), Copolyester (TPC-ES)

 – Medical device

 – DataMatrix Unique Device Identifier UDI

Professionals²

General instructions

Check the optimal and individual fit of the SecuTec Lumbo after the first fitting. Adapt the bridging frame (CFRP frame) according to the area of application (indications) and the build (anatomy) of the patient. Practice the correct procedure for donning of the product with the patient.

No improper modifications may be made to the product. Failure to comply may adversely affect the product's performance, thus canceling any product liability.

The SecuTec Lumbo is intended for the treatment of one patient only.

Part of the body this product is used for

In accordance with the indications (lower back). Please see intended purpose.

Assembly and fitting instructions

A – Knitted fabric portion
B – Fastening with finger pocket
C – Tensioning straps
D – Velour pocket for bridging frame (CFR frame)
E – Abdominal pad
F – Bridging frame (CFR frame)
G – Double-sided micro Velcro fastenings

Caution³

Wear protective gloves to avoid burns when handling the hot-air device.

To prevent damage to the material, please pay attention to the following: Use round bending hooks only. Do not heat the CFR frame at temperatures higher than 190°C. The optimum temperature for manipulating the frame is between 170 – 190°C. Maintain a minimum distance between the hot-air device and the CFRP frame. Heat the CFRP frame evenly and over the entire surface, not just at certain points.

Bend the CFRP frame in large arches only (no sharp bends). You can immerse the CFR frame in cold water to speed up the cooling process. This will not have any negative impact on the properties of the material.

You can heat and reshape the material more than once. The heating process does not visually alter the surface of the material. In order to find out the right time to shape the material, test it during the heating process.

❶ Heat the CFR frame with a hot-air device and shape it to produce the necessary bridging angle in accordance with the area of application (indication) and to suit the patient's anatomy. Make sure that the “up” marking on the frame is correctly oriented. Check the fit of the orthosis on the patient in both a standing and a sitting position after each modification.

❷ Lay out the velour pocket in front of you. The sewn-in label indicates which side should face upward “up”. Then turn the velour pocket so that the opening is pointing toward you. **❸** Place the CFR frame in the pocket and close it. **❹** Place the knitted section in front of you with the sewn-in label on the top right. Secure each pair of double-sided micro Velcro fastenings onto the velour strip on the knitted section. The fastenings must be spaced out at distances that enable the velour pocket to be securely connected to the frame and the knitted section. **❺** Fasten the velour pocket to the bridging frame and the knitted section.

❻ Fasten the abdominal pad onto the right-hand half of the fastening in front of you.

Practice the process of donning and doffing the orthosis with the patient.

Date of information: 2025-03

¹ Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso.

² A professional is any person who is authorized according to the state or institutional regulations as well as other provisions relevant to them for fitting and instruction in the use of orthopedic aids. As the manufacturer, we strongly recommend that this person be trained in orthopedics or possess comparable acquired competencies as well as dispose of orthopedic infrastructure.

³ Note on the risk of personal injury (risk of injury, harm to health, and accidents) or damage to property (damage to the product).

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