AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

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 english

Dear patient,

Thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please closely observe the specifications in these instructions for use and the instructions given by the professional².

If you have any questions, please contact your physician or medical retailer.

Intended purpose

SecuTec OA is a medical device. It is an orthosis¹ for relief and stabilization for osteoarthritis of the knee with leg axis misalignment.

Indications

- Medial or lateral osteoarthritis of the knee / arthritis (moderate to severe)
- Condition after cartilage transplant
- Following meniscus refixation
- Ruptured anterior / posterior cruciate ligament (ACL / PCL) with leg axis misalignment

Risks of using this product

 **Important information**

- The SecuTec OA orthosis is a prescribed product that must be worn under a physician's guidance. The product must be worn in accordance with these instructions for use and for the listed indications.
- No product liability is accepted if the product is not used or adjusted correctly.
- No improper modifications may be made to the product. Failure to comply with this requirement may adversely affect product performance, thereby excluding any product liability.
- Combining the orthosis with other products (e.g. compression stockings) must be discussed beforehand with your physician.
- Do not allow the product to come into contact with ointments, lotions, or substances containing grease or acids.
- After use in salt water or chlorinated water, please rinse the orthosis with plain water.
- Should any problems with the fit of the orthosis arise, for instance if the shape of your leg changes (e.g. due to an increase or reduction in muscle), or should the product have a functional defect, please contact a professional².
- Side effects involving the body as a whole have not been reported to date. Correct usage / donning of the product is assumed. Any aids – supports and orthoses – applied externally to the body can, if tightened excessively, lead to local pressure symptoms or, rarely, constrict blood vessels or nerves.
- For this reason, check the fit of the aid if you feel unpleasant pressure and consult the professional in this regard or for general questions about handling.
- If you notice any negative changes or an increase in symptoms while using the product, stop any further use and contact your physician immediately.

- Caution: when loosening the straps or tensioning and fastening systems, or when doffing the product or adjusting configured limitations of movement, there is not sufficient support / sufficient protection of the body part requiring treatment. In this case, particularly avoid placing any strain on this body part.
- Before using the product, especially following trauma or surgery on the knee joint or for immobilizing the leg, the professional or your physician must determine whether thrombosis prevention is necessary.
- Please adjust your behavior to the situation associated with your specific condition and any associated restriction of movement, and avoid anything that could delay recovery where possible. Proceed with particular caution.
- Furthermore, proceed with extreme caution, as your ability to stand and walk may be impacted.
- If you would like to use your Bauerfeind product due to acute symptoms / injuries, make sure to obtain and follow the advice of a professional before first using the product. Restrictions may exist given the specific indication which the professional can clarify. In particular, the ability to operate vehicles, other means of transport or machinery may be limited. In case of doubt, we recommend temporarily refraining from the activities mentioned above.
- If your perceptive faculties are restricted, if you take painkillers or blood-thinning medications or suffer from dementia, our products may only be worn when prescribed by a physician.

Contraindications

No clinically significant hypersensitivity reactions have been reported to date. If you have any of the following conditions, the aid should only be donned and worn after consultation with your physician:

1. Skin disorders / injuries in the relevant part of the body, particularly if inflammation is present. In case of any raised scars with swelling, redness or excessive heat build-up.
2. Varicose veins.
3. Impaired sensation and circulatory disorders of the legs / feet (e.g. diabetes mellitus).
4. Impaired lymphatic drainage – including unspecified soft tissue swelling.

Application instructions

Size calculation

The orthosis is selected and fitted by a professional after measuring the patient. You can refer to the size chart on the packaging for the correct size.

Procedure for the first fitting of the SecuTec OA knee orthosis

Adjustment options for the hinge (only by a professional²)

The hinge cover can be opened via a screw fastening. This cover fastens the movement-limiting wedges (extension – front / flexion – rear) and also protects the hinge mechanism against dirt.

Wedges are pre-installed at the factory. They primarily provide fixed stop absorption during extension and they do not limit the range of motion.

Opening the screw fastening (quarter turn until the arrow points at the "open" symbol  allows the hinge protection to be removed and the desired limitation blocks to be changed (➊ a – c).

• The desired extension and flexion wedges are selected in pairs from the optional frame and are inserted into the orthosis hinges on both sides. The following limitations are possible: Extension: 0°, 10°, 20°, 30°, 45° Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90° and fixation: 0°, 10°, 20°, 30°, 45°.

• Make sure that the screw fastening is showing the  symbol and click the cover back on to the hinge. When the wedges are flush with the surface of the cover, the screw fastening can be turned back to .

Adjustment (only by a professional²):

Preparation:

- The orthosis has special gel cushions on one side in order to ensure optimal pressure distribution.
- At the factory, the gel cushions for **correcting a varus knee (bandy legs)** or **medial osteoarthritis of the knee** are attached to the lateral hinge of the orthosis.
- In order to **correct a valgus knee (knock knees)** or **lateral osteoarthritis of the knee**, the gel cushions need to be exchanged.
- The oblong gel pad on the lateral side of the orthosis frame below the knee has a cut-out for the head of the fibula. Make sure that this is in the correct position.
- Mark the compromise pivot axis of the knee joint (as defined by Nietert) on the leg.
- The pivot point of the orthosis hinge lies in the center of the screw fastening on the hinge cover.
- When fitting the orthosis for the first time, start by unclicking all four quick-fit fasteners and loosening the two other straps.
- With the knee in a slightly flexed, unloaded position, don the SecuTec OA in such a way that the pivot point of the orthosis hinge corresponds to the compromise pivot axis (as defined by Nietert).

- Check the fit of the shell parts when they are laid flat to make sure their circumference and longitudinal direction are correct and reshape them if necessary. A bending iron with protective caps should be used for any extensive bending work.
- The calf shell is fitted to the anatomy of the lower leg by adjusting the slant.
- The position of the strap cushions can be altered as required.
- Now close the quick-fit fasteners that are close to the hinge (#1 and #2), followed by the ones further away (#3 and #4) ❶, and tighten the straps using the Velcro fastening ❷.
- Check once again that the orthosis hinges are in the correct position and tighten the Velcro straps close to the hinge (without a quick-fit fastener).
- Adjust the hinges on the lower leg frame to the shape of the leg; a slight correction position may be set here if required ❸.
- A 3 mm Allen key can be used for a valgus / varus setting of up to 25° in each case. The direction of rotation is marked on the orthosis ❹. Please note: extreme settings cause the patient pain and often offer no therapeutic benefit.
- Turn the screw on the hinge in the direction in which the lower part of the frame should move. A scale and a groove in the screw will help to guide you here. The distance between two marks corresponds to a correction of two degrees.
- Make the same setting on the other side in order to ensure that the shift is parallel; this results in an optimal application of force and ensures that the hinges are parallel ❺ (varus) ❻ (valgus).
- The gel condylar cushion can also be adjusted to the anatomy of the knee by using insertable wedges to increase the size of the contact surface for the correction. To make this adjustment, insert the wedges supplied (3 and/or 6 mm) from the front until the pad is lying flat against the knee ❼.
- **Please note:** in case of increased slipping of the orthosis despite adjustments, two additional siliconized strap pads can be ordered for replacement distally, close to the knee.

Donning of the SecuTec OA (for the patient)

- Since the straps and their length have been adjusted by the professional, these do not necessarily need to be undone. To don the orthosis, just undo the four quick-fit fasteners.
- Don the orthosis by holding it in front of you, placing your foot down through the frame, and then pulling the frame up to your knee ❶.
- Bend your leg slightly (approx. 30°). With your fingers, feel for the position of the kneecap. The center of the orthosis hinges is positioned at approximately the same height as the center of the kneecap.
- The grips of the quick-fit fasteners also have numbers (1–4) imprinted underneath indicating the best order for closing the straps ❷.
- Click all of the quick-fit fasteners on the straps into the eyelets, starting at the front with the fastener on the lower leg close to the hinge (#1), followed by the ones on the thigh close to the hinge (#2). Next, close the fasteners further away from the hinge, starting with the one on the lower leg first (#3), followed by the one on the thigh (#4).
- If the straps are too loose or too tight, you can tighten or loosen them slightly at this point.
- Attention: straps that are too loose may cause the orthosis to slip.
- After all the straps have been closed, the orthosis must lie evenly when the muscles are tensed.
- Once the orthosis has been put on, check once more that it is in the correct position ❸.

Doffing of the SecuTec OA

- Start by taking the quick-fit fasteners out of the eyelets on the back of the thigh frame of the orthosis.
- Next, undo the quick-fit fasteners of the straps at the front on the lower leg part.
- The front top strap and the rear bottom strap can remain closed.
- Once you have undone all of the quick-fit fasteners, push the orthosis down toward your foot and step out of the frame.

Cleaning instructions

Never expose the orthosis to direct heat (e.g. heater, sunlight, etc.)!

The aluminum parts of the orthosis are coated and can be cleaned using pH-neutral soap. We recommend washing the cushions and straps by hand at 30 °C. This product has been tested by our own quality control system. However, should you have any complaints, please contact your local medical supply retailer or the distributor from whom you purchased our product.

Part of the body this product is used for

In accordance with the indications (knee). Please see intended purpose.

Maintenance instructions

If handled and cared for correctly, the product is practically maintenance-free.

Assembly and fitting instructions

The SecuTec OA is supplied as a pre-assembled orthosis in standard sizes. The first donning, which includes the hinge adjustments and the static fitting, must be carried out by a professional.

Technical specifications / parameters

The SecuTec OA is an orthosis for the knee joint. It consists of a frame, hinges, straps, quick-fit fasteners and pads.

Notes on reuse

This product is intended to be fitted to one patient only.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. The warranty may be limited or excluded if the instructions on how to handle and care for the SecuTec OA have not been observed.

Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the professional were not observed
- The product was modified arbitrarily

Disclaimer

Do not diagnose yourself or decide on which medication to take unless you are a professional². Before first using our medical device, please seek the advice of a physician or professional² as this is the only way to assess the effect of our product on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of this professional as well as all instructions in this document or its online version, including any excerpts (as well as texts, images, graphics etc.). If you still have doubts after your consultation with the professional, please contact your physician or distributor, or contact us directly.

Duty to report

Due to regional legal regulations, you are obliged to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Aluminum (Al), Polyamide (PA), Zinc (Zn), Polyurethane (PUR), Styrene-ethylene/butylene-styrene block copolymers (TPS-SEBS), Polyoxymethylene (POM), stainless steel, Acrylate (AK), Polyethylene terephthalate (PET), Brass (CuZn), Polyester (PES), cork, steel

 – Medical device

 – DataMatrix Unique Device Identifier UDI

Date of information: 2025-03

¹ Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso.

² A professional is any person who is authorized according to the state or institutional regulations as well as other provisions relevant to them for fitting and instruction in the use of orthopedic aids. As the manufacturer, we strongly recommend that this person be trained in orthopedics or possess comparable acquired competencies as well as dispose of orthopedic infrastructure.

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