
Rev.1-2021-07_110408

de	deutsch	2	sk	slovensky	29
en	english	4	hu	magyar	31
fr	français	6	hr	hrvatski	33
nl	nederlands	9	sr	srpski	35
it	italiano	11	sl	slovenski	37
es	español	13	ro	romanian	40
pt	português	16	el	ελληνικά	42
no	norsk	18	tr	türkçe	44
fi	suomi	20	ru	русский	46
da	dansk	22	lv	latviešu	49
pl	polski	25	lt	lietuvių	51
cs	česky	27	uk	українська	53

de deutsch

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en english

Dear Customer,

thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please read through these **instructions for use** carefully. If you have any questions, please contact your doctor or medical retailer.

Intended purpose

SecuTec Omo is a medical device. It is an orthosis for immobilization of the shoulder with individual arm positioning and mobilization function.

The construction of the pelvic frame means that the orthosis can be adjusted and fitted for you. The design of the orthosis (fastening system, strap system) enables you to put the product on yourself.

The adjustable joint with several settings enables the arm to be positioned to suit the indications (abduction and rotation) and is therefore perfectly suited to everyday use. As a result, the SecuTec Omo has achieved a high level of compliance.

Indications

- Rotator cuff rupture (complex)
- Subcapital humeral head fracture
- SLAP lesion
- Shoulder replacement surgery
- Dislocated shoulder
- Post-traumatic, post-operative irritation
- Pre-operative care
- Acromioclavicular joint injury (AC rupture)

Risks of using this product

- ⚠ Caution***
- SecuTec Omo must only be worn in accordance with the specifications contained in these instructions for use and used for the areas of application (intended use, parts of the body) listed.
 - To ensure the best possible fit, the SecuTec Omo must be individually adjusted by a trained specialist¹. Only a trained specialist may carry out the (initial) fitting of the orthosis and provide instruction.
 - You should not alter the SecuTec Omo yourself. In the event of a product defect or wear to individual parts, please contact your medical retailer.
 - If you experience any unusual changes connected to the use of the product – e.g. an increase in discomfort – stop using it and contact your physician without delay.
 - The product may be damaged if it is not used correctly or is altered arbitrarily and, the warranty will not be valid.
 - Avoid contact with ointments, lotions, or substances containing grease or acids and keep the product away from direct sunlight.
 - Any aids applied externally to the body may, if tightened excessively, lead to local pressure points or, rarely, constrict the underlying blood vessels or nerves.
 - The use and duration of use must be discussed with the physician who prescribed the product. They will also inform you whether the orthosis can be used simultaneously with other products, and if so, which ones.

- Please note all the information you received with the product as well as the instructions for use.
- The sewn-in label of the SecuTec Omo, which contains information on the size, manufacturer, washing instructions and CE mark, is located on the inside of the abdominal strap.

Contraindications

- Side effects involving the body as a whole have not been reported to date. If you have any of the following conditions, the product should only be used after consultation with your physician:
- Skin disorders / injuries to the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up
 - Impaired sensation and circulatory disorders of the arms (e.g. diabetes mellitus)
 - Severe restriction of the performance of the heart and lungs (risk of a rise in blood pressure if the aid is worn during intense physical exertion)
 - Impaired lymph drainage – including soft tissue swellings of uncertain origin located away from the fitted aid
 - Mental and physical disabilities that render it impossible to use the orthosis safely

Application instructions

- A – Arm and hand support with Velcro straps
- B – Joint
- C – Pelvic frame with waist panels
- D – Abdominal fastening (panel) with finger pocket
- E – Velcro straps to adjust the circumference
- F – Abdominal strap
- G – Velcro tab
- H – Back fastening (panel)

Putting on the orthosis

The SecuTec Omo must be adapted by a trained specialist¹ before the first fitting. Wear the SecuTec Omo over a crease-free T-shirt or undershirt.

Ask someone to help you put on the orthosis if necessary.

1 Undo the abdominal strap. Unhook the abdominal fastening panel (with B) from the waist panel. Lift the round end slightly and pull it towards the arm to be treated. Undo the three straps on the arm support.

Once it is undone, position the SecuTec Omo on a chair without armrests. The arm support should be on the same side as the affected arm. The abdominal strap should be hanging down.

2 Sit between the two waist panels of the pelvic frame. Position the orthosis so that the recess in the pelvic frame sits directly above the iliac crest on the same side as the affected arm. **3** Put your arm into the arm support. **4** Fasten the abdominal strap. Hook the abdominal fastening panel into the waist panel. Pull the Velcro straps on the abdominal fastening panel tight. **5** Fasten all of the straps on the arm support. **6** Adjust the position of the orthosis when standing.

Putting on the arm support

7 Following consultation with your physician the arm support can be put on facing inward if necessary. To do this, press and hold down the blue adjustment button on the joint of the orthosis and position the affected arm with the arm support facing inward.

After releasing the adjustment button make sure that you hear the arm support engaging.

Using the SecuTec Omo with the shoulder strap

For an additional stabilizing effect, the optional shoulder strap can be used if necessary.

8 Attach the shoulder strap firmly to the outside of the SecuTec Omo on the back. Put the orthosis on following the instructions above. Guide the strap over the healthy shoulder to the front. Fasten the end of the shoulder strap with the square Velcro tab firmly to the abdominal strap. **9** Now undo the round Velcro end and adjust the length of the strap to fit.

Taking off the support

To take off the orthosis, start by undoing the abdominal fastening panel. Then undo all of the Velcro straps on the arm support.

Cleaning instructions

Note: Never expose the SecuTec Omo to direct heat (e.g. heater, sunlight, in your car, etc.!) This may damage the material, which can impair the effectiveness of the SecuTec Omo.

If necessary, you can wipe the plastic components of the SecuTec Omo with a soft cloth using a mild cleaning fluid. You can wash the pad by hand with a mild detergent. Remove the pad from the orthosis before washing. Leave the pad to air dry. **10** Reattach the pad using the snap fasteners.

Please refer to the care instructions on the sewn-in label of your SecuTec Omo. This can be found on the inside of the abdominal strap.

If you have any complaints, please contact your medical retailer only. Please note that only clean goods can be dealt with.

Part of the body this product is used for

In accordance with the indications (shoulder). See intended use.

Maintenance information

If handled and cared for correctly, the product is practically maintenance-free.

Assembly and fitting instructions

The SecuTec Omo is supplied as a pre-assembled orthosis in standard sizes and right and left versions. The position of the arm shell must be individually adjusted by a trained specialist.

Technical specifications / parameters, accessories

The SecuTec Omo is an orthosis for the shoulder. It is made up of an arm and hand support, a joint, a pelvic frame with two waist panels, and an abdominal strap with abdominal and back fastening panels.

Notes on reuse

This product is intended to be fitted to one patient only.

Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the SecuTec Omo have not been observed. Warranty is excluded if:

- The product was not used according to the indication
- The instructions given by the specialist were not observed
- The product was modified arbitrarily

Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

Disposal

Once you have finished using the product, please dispose of it according to local specifications.

Material content

Polyethylene, medium density (PE-MD), Polyamide glass fiber reinforced (PA-GF), Polyamide (PA), High Grade Steel, Polyurethane (PUR), Thermoplastic polyurethane (TPU), Polyester (PES), Steel, galvanized (Fe (Zn)), Brass, Polyoxymethylene (POM), Polypropylene (PP)

MD – Medical Device

UDI – DataMatrix Unique Device Identifier UDI

Specialists¹

General instructions

Adapt the orthosis according to the areas of application (indications) and the anatomy of the patient. Check the optimal and individual fit of the SecuTec Omo after the first fitting. Practice the correct procedure for putting on the orthosis with the patient. The product must not be changed in any improper way. Failure to comply with this requirement may adversely affect product performance, thereby voiding any product liability. The SecuTec Omo is intended for the treatment of one patient only.

Part of the body this product is used for

In accordance with the indications (shoulder).
See intended use.

Assembly and fitting instructions

- A – Arm and hand support with Velcro straps
- B – Joint
- C – Pelvic frame with waist panels
- D – Abdominal fastening (panel) with finger pocket
- E – Velcro straps to adjust the circumference
- F – Abdominal strap
- G – Velcro tab
- H – Back fastening (panel)

Tools:

Included in the package:
Screw extension (supplied attached to the lower wing screw)
Not included in the package: Scissors (for plastic)

1 Undo the abdominal strap. Unhook the abdominal fastening panel (with B) from the waist panel. Lift the round end slightly and pull it towards the arm to be treated. Separate the abdominal fastening panel from the strap at the wide Velcro tab. Now undo the three straps on the arm support.

2 Position the orthosis so that the recess in the pelvic frame sits directly above the iliac crest on the same side as the affected arm.

3 Place the arm on the arm and hand support. Make sure that the joint for setting the rotation is vertically level with the shoulder joint. Shortening the waist panels

4 Now check the length of the waist panels. If the waist measurement is less than 80 cm, shorten the waist panels of the pelvic shell so that the distance between the ends of the front and back waist panels is at least 15 cm after the orthosis is put on. It is better to shorten the front waist panel. (If the back waist panel is shortened, the edges may press into the back when lying down.) To do this, first cut the plastic along the markings, then shorten the pad, leaving at least 1 cm protruding over the edge of the plastic. If necessary, the pad can also be removed completely for shortening by undoing the fastening clips.

5 Now guide the abdominal strap behind the back to the opposite side of the pelvis.

6 To ensure that the waist panels and abdominal strap of the orthosis will always be positioned horizontally on the patient when standing, regardless of the shape of the pelvis, the slope of the waist panels can be adjusted. To do this, loosen the unions on both waist panels, adjust the angle until they are horizontal and tighten the unions again afterward. Use the enclosed Allen key for this.

Shortening the abdominal strap

7 Shorten the abdominal strap, leaving a distance of approx. 20 cm between the strap and the front waist panel after the strap has been cut. The strap can also be cut diagonally to fit the anatomy (pelvis shape) better.

8 Attach the abdominal fastening panel to the shortened abdominal strap.

9 Fasten the abdominal strap by hooking the fastening panel to the waist panel. Pull both adjustment straps tight. Setting the abduction angle and adjusting the arm support to the length of the arm

10 Loosen the upper locking screw, use it to set the abduction angle determined by the physician and tighten the screw again. Attach the covering caps to the front and back of the arm support.

11 Now fit the arm support to the arm. Loosen the lower locking screw using the attached screw extension. Move the arm support up so that the arm is lying flush against the support and the elbow is in the cut-out provided. When doing this, make sure that both shoulders are at the same height. Then tilt the arm support. To avoid neurological impairment or deficits, the edges of the arm support must not put pressure on the arm.

12 Secure the setting by tightening the lower locking screw. Pull off the screw extension. Attach the covering caps to the front and back of the arm support.

13 Now fasten the Velcro straps on the arm support and cut off the excess with scissors. When doing this, make sure that thicker clothing could still be worn under the orthosis.

Adjusting the hand support

14 Adjust the hand support to the length of the forearm by pushing it forward or backward until the palm lies flush against the support.

If necessary, the hand support can be removed completely. Setting the outward and inward rotation on the arm support Following consultation with the physician, the patient can adjust the outward and inward rotation of the arm in the arm support. To enable this, release or lock the required direction of movement.

15 Outward rotation is locked by default; it can be released if necessary by removing the locking clip. Inward rotation is enabled by default; it can be locked if necessary by inserting the locking clip.

16 Practice putting on the arm support with the patient: The patient should hold down the blue adjustment button on the joint of the orthosis and position the affected arm with the arm support facing inward.

Since the product is usually fitted directly after surgery and worn around the clock for several weeks, the settings should be checked again after a few hours or after the first night and readjusted.

Shoulder strap

17 Attach the straight end of the shoulder strap to the fastening panel on the back. Guide the strap diagonally over the healthy shoulder. Attach the adjustment clasp to the abdominal fastening panel and thread the shoulder strap through the loop. Pull the shoulder strap tight and attach it to the abdominal strap.

Version: 2021-07

* Note on the risk of personal injury (risk of injury, harm to health and accidents) or damage to property (damage to the product).

1 A specialist is any person who is authorized according to the state regulations for fitting and instruction in the use of orthoses which are relevant to you.

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